


PERISTALTIC PUMP

6-channel Pump



Features

- Can deliver up to six kinds of liquids simultaneously, despite a compact size.
- Easy replacement of tubes
- Stepping motor type also available, which can handle microliter flow rates.
- Observed flow rate accuracy is $\pm 10\%$ using six tubes from the same lot.

Specifications

Motor	DC Motor	Bipolar Stepper Motor
Flow Rate	830 $\mu\text{l}/\text{min}$	0.23 - 350 $\mu\text{l}/\text{min}$
Tubing Material	Silicone	
Tubing Diameter	I.D. 1.0 \times O.D. 2.0 mm	
Pump Pressure	30 kPa	
Rated Voltage	DC 3 V	DC 10 V
Dimensions	31 x 84 x 32 mm	

Note: Specifications etc. may be changed at any time without notice.

TFS (Takasago Fluidic Systems)

1900 West Park Drive, Suite 280, Westborough, MA 01581 U.S.A.

Tel: 508-983-1434 FAX: 508-983-1401

E-mail: info@takasago-fluidics.com URL: <http://www.takasago-fluidics.com/>